INTENDED USE

QConnect Blue is a multi-marker quality control (QC) sample intended to monitor the performance of in-vitro assays for determination of specific analytes in human serum and plasma, as listed in "Contents" below. QConnect Blue samples are optimized for assays and platforms including ABBOTT ARCHITECT, ABBOTT Alinity s and DiaSorin LIAISON and are unassayed controls without target values. These controls must not be used as a substitute for the mandatory manufacturer's kit controls provided with the assay. Examples of reactivity (no target values) for this control on different test systems can be found at <http://www.diamex.com>.

QConnect QC results are helpful in determining the precision of testing systems and in identifying sources of variation. Users of QConnect QC samples have access to EDCNet, an Internet-based QC results monitoring system at <http://www.nrlquality.org.au/qconnect>.

CONTENTS

QConnect Blue contains human plasma with sodium azide and ProClin 950 as preservative. QConnect Blue will provide a positive reaction in assays that detect the following analytes:

Anti-HIV 1 IgG	Anti-HBc	Anti-HTLV IgG*
Anti-HCV IgG	HBsAg	Anti- <i>Treponema pallidum</i> IgG*

*Detection of this analyte is currently not intended for ABBOTT Alinity s.

The negative control QConnect SeroNeg will provide negative/non-reactive results to the above mentioned analytes. QConnect Blue was confirmed negative for HBV DNA, HIV RNA and HCV RNA using a PCR method.

STORAGE AND STABILITY

QConnect Blue:

- Can be used until the expiration date when stored at 2-8 °C;
- Is stable for 90 days once opened providing it is closed tightly after use and returned to storage at 2-8 °C;
- Should always be stored upright;
- Must not be used beyond the expiration date.

INSTRUCTIONS FOR USE

- QConnect Blue is ready-to-use;
- Mix thoroughly by gentle swirling prior to use;
- For automated analyzers, label the vial with the barcode provided and then place the QConnect QC in a sample rack as per routine samples;
- Return vials to 2-8 °C storage after use.

LIMITATIONS

- Do not use this product beyond the expiration date;
- Do not use this product where there is evidence of microbial contamination and/or high turbidity;
- Use QConnect Blue as per assay manufacturer's instructions for use only;
- Deviation from the recommended procedure may lead to unreliable results;
- This product is provided for quality assurance only and should not be used for calibration purposes;
- Possible causes of variation in QConnect QC performance:
 - Different assay reagent lots;
 - Different calibrator lots used;
 - Different analyzers used;
 - Different operators;
 - Time of day when testing.

PRECAUTIONS AND DISPOSAL INFORMATION

- QConnect Blue is for in-vitro diagnostic use only and should be tested by trained personnel only;
- QConnect Blue contains materials of human origin and should be considered potentially infectious. Observe universal precautions for prevention of transmission of infectious agents when handling this product;
- Do not pipette by mouth;
- The preservative sodium azide may react with copper or lead to form dangerous compounds. Although QConnect Blue contains only a small amount of sodium azide, please dilute with a sufficient amount of water when disposing these products into waste water systems;
- This product also contains methylisothiazolones, which are a compound of ProClin. Methylisothiazolones are classified per applicable European Community (EC) Directives as: skin sensitization;
- Discard all material in a safe manner that is in compliance with local and national regulations and in accordance to laboratory procedure for discard of clinical waste.

! WARNING: Contains methylisothiazolones.
H317 May cause an allergic skin reaction.

Prevention:

- P261 Avoid breathing mist/ vapours/ spray.
- P272 Contaminated work clothing should not be allowed out of the workplace.
- P280 Wear protective gloves/ protective clothing.

Response:

- P302 + P352 IF ON SKIN: Wash with plenty of soap and water.
- P333 + P313 If skin irritation or rash occurs: Get medical advice/ attention.
- P363 Wash contaminated clothing before reuse.

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