## INTENDED USE

QConnect Mumps M is a single marker quality control (QC) sample intended to monitor the performance of in-vitro assays for determination of anti-mumps virus IgM antibodies in human serum and plasma. QConnect Mumps M samples are optimized for assays and platforms including DiaSorin LIAISON and Siemens Enzygnost and is an unassayed control without target values. This control must not be used as a substitute for the mandatory manufacturer's kit controls provided with the assay. Examples of reactivity (no target values) for this control on different test systems can be found at <http://www.diamex.com>.

QConnect QC results are helpful in determining the precision of testing systems and in identifying sources of variation. Users of QConnect QC samples have access to EDCNet, an Internet-based QC results monitoring system at <http://www.nrlquality.org.au/qconnect>.

## CONTENTS

QConnect Mumps M contains human plasma with sodium azide and ProClin 950 as preservatives. Each plasma unit was separately tested and found negative/non-reactive for anti-HIV 1/2, anti-HCV and HBsAg.

QConnect Mumps M will provide a positive reaction with the following analyte:

Anti-mumps virus IgM

The negative control QConnect Pediatric M will provide negative/non reactive results to the above mentioned analyte.

## STORAGE AND STABILITY

QConnect Mumps M:

- Can be used until the expiration date when stored at 2-8 °C;
- Is stable for 90 days once opened providing it is closed tightly after use and returned to storage at 2-8 °C;
- Should always be stored upright;
- Must not be used beyond the expiration date.

## INSTRUCTIONS FOR USE

- QConnect Mumps M is ready-to-use;
- Mix thoroughly by gentle swirling prior to use;
- For automated analyzers, label the vial with the barcode provided and then place the QConnect QC in a sample rack as per routine samples;
- Return vials to 2-8 °C storage after use.

## LIMITATIONS

- Do not use this product beyond the expiration date;
- Do not use this product where there is evidence of microbial contamination and/or high turbidity;
- Use QConnect Mumps M as per assay manufacturer's instructions for use only;
- Deviation from the recommended procedure may lead to unreliable results;
- This product is provided for quality assurance only and should not be used for calibration purposes;
- Possible causes of variation in QConnect QC performance:
  - Different assay reagent lots;
  - Different calibrator lots used;
  - Different analyzers used;
  - Different operators;
  - Time of day when testing.

## PRECAUTIONS AND DISPOSAL INFORMATION

- QConnect Mumps M is for in-vitro diagnostic use only and should be tested by trained personnel only;
- QConnect Mumps M contains materials of human origin and should be considered potentially infectious. Observe universal precautions for prevention of transmission of infectious agents when handling this product;
- Do not pipette by mouth;
- The preservative sodium azide may react with copper or lead to form dangerous compounds. Although QConnect Mumps M contains only a small amount of sodium azide, please dilute with a sufficient amount of water when disposing these products into waste water systems;
- This product also contains methylisothiazolones, which are a compound of ProClin. Methylisothiazolones are classified per applicable European Community (EC) Directives as: skin sensitization;
- Discard all material in a safe manner that is in compliance with local and national regulations and in accordance to laboratory procedure for discard of clinical waste.

**WARNING:** Contains methylisothiazolones.  
H317 May cause an allergic skin reaction.

## Prevention:

- P261 Avoid breathing mist/ vapours/ spray.  
P272 Contaminated work clothing should not be allowed out of the workplace.  
P280 Wear protective gloves/ protective clothing.

## Response:

- P302 + P352 IF ON SKIN: Wash with plenty of soap and water.  
P333 + P313 If skin irritation or rash occurs: Get medical advice/ attention.  
P363 Wash contaminated clothing before reuse.

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GmbH
DiaMex GmbH, Siemensstr. 38, 69123 Heidelberg, Germany, BE134ML, 2018-05-23