**INTENDED USE**

QConnect ToRCH G is a multi-marker quality control (QC) sample intended to monitor the performance of in-vitro assays for determination of specific analytes in human serum and plasma, as listed in "Contents" below. QConnect ToRCH G samples are optimized for assays and platforms including ABBOTT ARCHITECT, DiaSorin LIAISON and Siemens Enzygnost and are unassayed controls without target values. These controls must not be used as a substitute for the mandatory manufacturer's kit controls provided with the assay. Examples of reactivity (no target values) for this control on different test systems can be found at <http://www.diamex.com>.

QConnect QC results are helpful in determining the precision of testing systems and in identifying sources of variation. Users of QConnect QC samples have access to EDCNet, an Internet-based QC results monitoring system at <http://www.nrlquality.org.au/qconnect>.

**CONTENTS**

QConnect ToRCH G contains human plasma with sodium azide and ProClin 950 as preservative. Each plasma unit was separately tested and found negative/non-reactive for anti-HIV1/2, anti-HCV and HBsAg. The positive control contains IgG-antibodies to:

<i>Toxoplasma gondii</i>	Rubella-virus	Cytomegalovirus
Epstein-Barr-virus VCA	Epstein-Barr-virus EBNA	Herpes-simplex-virus 1/2

The negative control QConnect ToRCH G will provide negative/non-reactive results to the above mentioned analytes.

**STORAGE AND STABILITY**

QConnect ToRCH G:

- Can be used until the expiration date when stored at 2-8 °C;
- Is stable for 90 days once opened providing it is closed tightly after use and returned to storage at 2-8 °C;
- Should always be stored upright;
- Must not be used beyond the expiration date.

**INSTRUCTIONS FOR USE**

- QConnect ToRCH G is ready-to-use;
- Mix thoroughly by gentle swirling prior to use;
- For automated analyzers, label the vial with the barcode provided and then place the QConnect QC in a sample rack as per routine samples;
- Return vials to 2-8 °C storage after use.

**LIMITATIONS**

- Do not use this product beyond the expiration date;
- Do not use this product where there is evidence of microbial contamination and/or high turbidity;
- Use QConnect ToRCH G as per assay manufacturer's instructions for use only;
- Deviation from the recommended procedure may lead to unreliable results;
- This product is provided for quality assurance only and should not be used for calibration purposes;
- Possible causes of variation in QConnect QC performance:
  - Different assay reagent lots;
  - Different calibrator lots used;
  - Different analyzers used;
  - Different operators;
  - Time of day when testing.

**PRECAUTIONS AND DISPOSAL INFORMATION**

- QConnect ToRCH G is for in-vitro diagnostic use only and should be tested by trained personnel only;
- QConnect ToRCH G contains materials of human origin and should be considered potentially infectious. Observe universal precautions for prevention of transmission of infectious agents when handling this product;
- Do not pipette by mouth;
- The preservative sodium azide may react with copper or lead to form dangerous compounds. Although QConnect ToRCH G contains only a small amount of sodium azide, please dilute with a sufficient amount of water when disposing these products into waste water systems;
- This product also contains methylisothiazolones, which are a compound of ProClin. Methylisothiazolones are classified per applicable European Community (EC) Directives as: skin sensitization;
- Discard all material in a safe manner that is in compliance with local and national regulations and in accordance to laboratory procedure for discard of clinical waste.

 **WARNING:** Contains methylisothiazolones. H317 May cause an allergic skin reaction.

**Prevention:**

- P261 Avoid breathing mist/ vapours/ spray.
- P272 Contaminated work clothing should not be allowed out of the workplace.
- P280 Wear protective gloves/ protective clothing.

**Response:**

- P302 + P352 IF ON SKIN: Wash with plenty of soap and water.
- P333 + P313 If skin irritation or rash occurs: Get medical advice/ attention.
- P363 Wash contaminated clothing before reuse.

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