



## Anti-SARS-CoV-2 Tests Manufacturer independent Optitrol COVID-19 IgG and IgM Controls

- Human plasma-based control treated the same as patient specimen
- Assessment of laboratory precision/accuracy of testing systems
- Third party independent control verifying systematic errors
- Satisfies regulatory guidelines for use of independent controls
- High lot to lot consistency
- Reacts at an appropriate level within the analytical range of the Anti-SARS-CoV-2 Tests
- Ready-to-use and bar-coded to ensure traceability
- Electronic Instructions For Use (eIFU)
- Track product performance with EDCNet
- Storage and Stability  
Sealed: 12 months at 2 to 8°C  
After opening: 60 days at 2 to 8°C



### Ordering information:

DiaMex GmbH  
Siemensstraße 38  
69123 Heidelberg  
[order@diamex.com](mailto:order@diamex.com)

Product code	QC Product	Kit Format	Vial Volume [mL]
SR14039	Optitrol Covid-19 IgG	1 x pos	1 mL
SR14037	Optitrol Covid-19 IgG	2 x pos.	2,5 mL
SR14033	Optitrol Covid-19 IgG	4 x pos.	2,5 mL
SR14049	Optitrol Covid-19 IgM	1 x pos	1 mL
SR14047	Optitrol Covid-19 IgM	2 x pos.	2,5 mL
SR14043	Optitrol Covid-19 IgM	4 x pos.	2,5 mL