

Representative Reactivities Table

Optitrol NAT HAV



Positive Control

Version: 1 Date: June 16, 2025

Analyte	Manufacturer	Assay Name	Unit	Mean	SD	Samples	Labs	Representative Reactivity
HAV	Roche	cobas® dpx for use on cobas® 6800 systems	Ct	30.41	0.15	10	1	positive
	Altona	AltoStar® HAV RT-PCR Kit 1.5	Ct	33.06	0.16	10	1	positive

The stated mean values were obtained by testing one product batch. They show possible representative values but no target values.
Test results may vary among used instruments, assays, reagent lots, and different laboratory protocols.