

Representative Reactivities Table



Optitrol Syphilis G

Positive Control

Date: April 26, 2018

Analyte	Manufacturer	Assay Name	Method	Unit	Mean	SD	Samples	Labs	Representative Reactivity
Anti-<i>Treponema pallidum</i> IgG	Abbott	Architect Syphilis TP	CMIA	S/CO	11.8	0.25	10	2	positive
	Fujirebio	Serodia <i>Treponema pallidum</i> particle agglutination (TP-PA) test	PA	Titer	> 1:320	n.a.	7	1	positive
	Newmarket Biomedical	Newbio-PK TPHA	HA	SPC	3.5	n.a.	12	2	positive
Anti-"reagin" antibody	Biokit	RPR retest	Rapid test	Titer	1:2	n.a.	7	1	positive

n.a.: not applicable

The stated mean values were obtained by testing one product batch. They show possible representative values but no target values. Test results may vary among used instruments, assays, reagent lots, and different laboratory protocols.