

Representative Reactivities Table



Optitrol ToRCH M neg

Negative Control

Date: April 06, 2021

Analyte	Manufacturer	Assay Name	Method	Representative Reactivity
Rubella IgM	Abbott	Architect Rubella IgM	CMIA	negative
Toxo IgM	Abbott	Architect Toxo IgM	CMIA	negative
CMV IgM	Abbott	Architect CMV IgM	CMIA	negative
EBV VCA IgM	Abbott	Architect EBV VCA IgM	CMIA	negative
Rubella IgM	DiaSorin	Liaison Rubella IgM	CLIA	negative
Toxo IgM	DiaSorin	Liaison Toxo IgM	CLIA	negative
CMV IgM	DiaSorin	Liaison CMV IgM	CLIA	negative
EBV VCA IgM	DiaSorin	Liaison EBV/VCA IgM	CLIA	negative
HSV 1/2 IgM	DiaSorin	Liaison HSV 1/2 IgM	CLIA	negative

The stated mean values were obtained by testing one product batch. They show possible representative values but no target values.
Test results may vary among used instruments, assays, reagent lots, and different laboratory protocols.