

Representative Reactivities Table



Optitrol Paediatric G-2

Positive Control

Version: 1

Date: May 24, 2022

Analyte	Manufacturer	Assay Name	Method	Unit	Mean	SD	Samples	Labs	Representative Reactivity
Anti-Measles Virus IgG	bioMérieux	VIDAS Measles IgG	ELFA	Test value	2.7	0.14	10	1	positive
Anti-Mumps Virus IgG	bioMérieux	VIDAS Mumps IgG	ELFA	Test value	1.6	0.18	10	1	positive
Anti-Varicella Zoster Virus IgG	bioMérieux	VIDAS Varicella-Zoster IgG	ELFA	Test value	3.0	0.17	10	1	positive

The stated mean values were obtained by testing one product batch. They show possible representative values but no target values. Test results may vary among used instruments, assays, reagent lots, and different laboratory protocols.