

Representative Reactivities Table



Optitrol Paediatric G

Positive Control

Version: 3

Date: May 18, 2022

Analyte	Manufacturer	Assay Name	Method	Unit	Mean	SD	Samples	Labs	Representative Reactivity
Anti-Measles Virus IgG	DiaSorin	LIAISON Measles IgG	CLIA	AU/mL	61.6	1.9	10	2	positive
Anti-Mumps Virus IgG	DiaSorin	LIAISON Mumps IgG	CLIA	AU/mL	50.9	2.4	10	2	positive
Anti-Varicella Zoster Virus IgG	DiaSorin	LIAISON VZV IgG	CLIA	mIU/mL	366	11	10	2	positive
Anti-Parvovirus B19 IgG	DiaSorin	LIAISON Biotrin Parvovirus B19 IgG	CLIA	Index	4.50	0.19	10	2	positive

The stated mean values were obtained by testing one product batch. They show possible representative values but no target values. Test results may vary among used instruments, assays, reagent lots, and different laboratory protocols.

