

Representative Reactivities Table



Optitrol ToRCH G neg

Negative Control

Date: April 06, 2021

Analyte	Manufacturer	Assay Name	Method	Representative Reactivity
Rubella IgG	Abbott	Architect Rubella IgG	CMIA	negative
Toxo IgG	Abbott	Architect Toxo IgG	CMIA	negative
CMV IgG	Abbott	Architect CMV IgG	CMIA	negative
EBV VCA IgG	Abbott	Architect EBV VCA IgG	CMIA	negative
EBNA IgG	Abbott	Architect EBNA-1 IgG	CMIA	negative
Rubella IgG	DiaSorin	Liaison Rubella IgG	CLIA	negative
Toxo IgG	DiaSorin	Liaison Toxo IgG II	CLIA	negative
CMV IgG	DiaSorin	Liaison CMV IgG II	CLIA	negative
EBV VCA IgG	DiaSorin	Liaison VCA IgG	CLIA	negative
HSV-1/2 IgG	DiaSorin	Liaison HSV-1/2 IgG	CLIA	negative
EBNA IgG	DiaSorin	Liaison EBNA IgG	CLIA	negative
Rubella IgG	Siemens	Enzygnost Anti- Rubella/IgG	ELISA	negative
Toxo IgG	Siemens	Enzygnost Toxo/IgG	ELISA	negative
CMV IgG	Siemens	Enzygnost Anti-CMV/IgG	ELISA	negative
EBV VCA IgG	Siemens	Enzygnost Anti-EBV IgG	ELISA	negative
EBNA IgG	Siemens	Enzygnost Anti-EBNA/IgG	ELISA	negative
HSV-1/2 IgG	Siemens	Enzygnost HSV-1/2/IgG	ELISA	negative

The stated mean values were obtained by testing one product batch. They show possible representative values but no target values. Test results may vary among used instruments, assays, reagent lots, and different laboratory protocols.