

Representative Reactivities Table



QConnect Measles M

Positive Control

Version: 1

Date: May 25, 2018

Analyte	Manufacturer	Assay Name	Method	Unit	Mean	SD	Samples	Labs	Representative Reactivity
Anti-measles virus IgM	DiaSorin	LIAISON Measles IgM	CLIA	Index	2.28	0.18	11	2	positive
	Siemens	Enzygnost Anti-Measles Virus/IgM	ELISA	Index / Ratio	2.97	1.10	11	2	positive

The stated mean values were obtained by testing one product batch. They show possible representative values but no target values. Test results may vary among used instruments, assays, reagent lots, and different laboratory protocols.