

## Representative Reactivities Table



### QConnect Mumps M

Positive Control

Version: 1

Date: May 25, 2018

Analyte	Manufacturer	Assay Name	Method	Unit	Mean	SD	Samples	Labs	Representative Reactivity
Anti-mumps virus IgM	DiaSorin	LIAISON Mumps IgM	CLIA	Index	2.36	0.41	14	2	positive
	Siemens	Enzygnost Anti-Paroditis-Virus/IgM	ELISA	Index / Ratio	2.69	0.57	10	2	positive

The stated mean values were obtained by testing one product batch. They show possible representative values but no target values. Test results may vary among used instruments, assays, reagent lots, and different laboratory protocols.