

## Representative Reactivities Table



### QConnect Parvo M

Positive Control

Version: 1

Date: May 25, 2018

Analyte	Manufacturer	Assay Name	Method	Unit	Mean	SD	Samples	Labs	Representative Reactivity
Anti-parvovirus B19 IgM	DiaSorin	LIAISON Biotrin Parvovirus B19 IgM	CLIA	Index	4.85	0.23	11	2	positive
	NovaTec	Novagnost Parvovirus B19 IgM	ELISA	NTU	9.39	0.74	10	2	negative/grey zone

The stated mean values were obtained by testing one product batch. They show possible representative values but no target values. Test results may vary among used instruments, assays, reagent lots, and different laboratory protocols.